

ImmunoCAP™ Total IgE

Fluorenzymeimmunoassay

CLIA Complexity Category = Moderately Complex

Directions for Use 52-5292-EN/08

INTENDED USE

ImmunoCAP Total IgE is an in vitro test system for the quantitative measurement of circulating total IgE in human serum or plasma. It is intended for in vitro diagnostic use as an aid in the clinical diagnosis of IgE mediated allergic disorders in conjunction with other clinical findings, and is to be used in clinical laboratories. ImmunoCAP Total IgE is to be used with the instruments Phadia 100, Phadia 200, Phadia 250, Phadia 1000, Phadia 2500 or Phadia 5000.

SUMMARY AND EXPLANATION OF THE TEST

Since 1967, when the first assays for serum immunoglobulin E (I-3) were described, these measurements have become well established components of the investigation of allergic patients. The serum concentration of IgE is significantly elevated in patients suffering from extrinsic asthma, hayfever or atopic eczema. The increase during childhood is slow. Adult values are not stabilized until 15-20 years of age (6-7).

PRINCIPLE OF THE PROCEDURE

Anti-IgE, covalently coupled to ImmunoCAP, reacts with the total IgE in the patient sample. After washing, enzyme labeled antibodies against IgE are added to form a complex. Following incubation, unbound enzyme-anti-IgE is washed away and the bound complex is then incubated with a developing agent. After stopping the reaction, the fluorescence of the eluate is measured. The fluorescence is directly proportional to the concentration of IgE in the sample. The higher the response, the more IgE is present in the sample. To evaluate the test results, the responses for the patient samples are transformed to concentrations with the use of a calibration curve.

REAGENTS AND MATERIAL

Reagents are packaged as described below, each purchased separately. The two digit suffix (-xx) on the article number may vary between countries. All kits are not available in all countries. The expiration date and storage temperature are stated on the labels. Do not use reagents beyond their expiration dates.

Note: It is not recommended to pool any reagents.

Keep the ImmunoCAP carrier closed to avoid evaporation of buffer.

Reagents for Phadia 100

- **ImmunoCAP Total IgE** (Art No 10-9517-01: for 96 determinations)
 - Total IgE Conjugate (1 vial)
 - Total IgE Curve Control 1 (CC-1) (4 single dose vials)
 - Total IgE Curve Control 2 (CC-2) (4 single dose vials)
- **ImmunoCAP Total IgE Curve Controls** (CC-1 and CC-2) (Art No 10-9257-01: for 3 additional assay runs)
- **ImmunoCAP Total IgE Calibrators** (Cal-xx) (Art No 10-9252-01: for 1 calibration curve)
- **ImmunoCAP Total IgE Anti-IgE** (a-IgE) (Art No 14-4509-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9478-01: for 600 determinations)
- **Stop Solution** (Art No 10-9479-01: for 600 determinations)
- **Washing Solution** (Art No 10-9422-01: 6 x 1 l)
 - Washing Solution Additive, 6 x 17.2 ml
 - Washing Solution Concentrate, 6 x 80 ml
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml
- **ImmunoCAP Total IgE Control LMH** (Art No 10-9447-01: for 3 x 2 x 4 determinations)

Reagents for Phadia 200

- **ImmunoCAP Total IgE Conjugate 100** (Art No 10-9319-01: for 6 x 100 determinations)
- **ImmunoCAP Total IgE Conjugate 400** (Art No 10-9480-01: for 6 x 400 determinations)
- **ImmunoCAP Total IgE Curve Control Strip** (CC-1 and CC-2) (Art No 10-9325-01: for 5 x 3 curve controls)
- **ImmunoCAP Total IgE Calibrator Strip** (Art No 10-9387-01: for 5 calibration curves)
- **ImmunoCAP Total IgE Anti-IgE** (a-IgE) (Art No 14-4509-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9441-01: for 6 x 200 determinations; Art No 10-9440-01: for 6 x 315 determinations)
- **Stop Solution** (Art No 10-9442-01: for 6 x 185 determinations; Art No 10-9479-01: for 6 x 100 determinations)
- **Washing Solution** (Art No 10-9422-01: 6 x 1 l)
 - Washing Solution Additive, 6 x 17.2 ml
 - Washing Solution Concentrate, 6 x 80 ml
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml

- **ImmunoCAP Total IgE Control LMH** (Art No 10-9447-01: for 3 x 2 x 4 determinations)

Reagents for Phadia 250

- **ImmunoCAP Total IgE Conjugate 100** (Art No 10-9319-01: for 6 x 100 determinations)
- **ImmunoCAP Total IgE Conjugate 400** (Art No 10-9480-01: for 6 x 400 determinations)
- **ImmunoCAP Total IgE Curve Control Strip** (CC-1 and CC-2) (Art No 10-9325-01: for 5 x 3 curve controls)
- **ImmunoCAP Total IgE Calibrator Strip** (Art No 10-9387-01: for 5 calibration curves)
- **ImmunoCAP Total IgE Anti-IgE** (a-IgE) (Art No 14-4509-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9441-01: for 6 x 200 determinations; Art No 10-9440-01: for 6 x 315 determinations)
- **Stop Solution** (Art No 10-9442-01: for 6 x 185 determinations)
- **Washing Solution** (Art No 10-9422-01: 6 x 1 l)
 - Washing Solution Additive, 6 x 17.2 ml
 - Washing Solution Concentrate, 6 x 80 ml
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml
- **ImmunoCAP Total IgE Control LMH** (Art No 10-9447-01: for 3 x 2 x 4 determinations)

Reagents for Phadia 1000

- **ImmunoCAP Total IgE Conjugate 100** (Art No 10-9319-01: for 6 x 100 determinations)
- **ImmunoCAP Total IgE Conjugate 400** (Art No 10-9480-01: for 6 x 400 determinations)
- **ImmunoCAP Total IgE Curve Control Strip** (CC-1 and CC-2) (Art No 10-9325-01: for 5 x 3 curve controls)
- **ImmunoCAP Total IgE Calibrator Strip** (Art No 10-9387-01: for 5 calibration curves)
- **ImmunoCAP Total IgE Anti-IgE** (a-IgE) (Art No 14-4509-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9439-01: for 6 x 1200 determinations; Art No 10-9314-01: for 6 x 2000 determinations)
- **Stop Solution** (Art No 34-2271-51: for 1200 determinations)
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml
- **ImmunoCAP Total IgE Control LMH** (Art No 10-9447-01: for 3 x 2 x 4 determinations)

Reagents for Phadia 2500 and Phadia 5000

- **ImmunoCAP Total IgE Conjugate 100** (Art No 10-9319-01: for 6 x 100 determinations)
- **ImmunoCAP Total IgE Conjugate 400** (Art No 10-9480-01: for 6 x 400 determinations)
- **ImmunoCAP Total IgE Curve Control Strip** (CC-1 and CC-2) (Art No 10-9325-01: for 5 x 3 curve controls)
- **ImmunoCAP Total IgE Calibrator Strip** (Art No 10-9387-01: for 5 calibration curves)
- **ImmunoCAP Total IgE Anti-IgE** (a-IgE) (Art No 14-4509-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9314-01: for 6 x 2000 determinations)
- **Stop Solution** (Art No 34-2337-11: for 4600 determinations)
- **Washing Solution Additive** (Art No 10-9518-01: 4 x 850 ml)
- **Washing Solution Concentrate** (Art No 34-2337-21: 1 x 2800 ml)
- **ImmunoCAP Total IgE Control LMH** (Art No 10-9447-01: for 3 x 2 x 4 determinations)

Details of reagents

ImmunoCAP Total IgE Conjugate/Conjugate 100/Conjugate 400	
β-Galactosidase-anti-IgE Approximately 2 µg/ml (mouse monoclonal antibodies) Sodium azide 0.06%	Ready for use. Store at 2 – 8 °C until expiration date. Do not freeze!

ImmunoCAP Total IgE Curve Controls/Curve Control Strip	
(human IgE in buffer) Preservative* <0.003%	Ready for use. Store at 2 – 8 °C until expiration date.

ImmunoCAP Total IgE Calibrators/Calibrator Strip	
(human IgE in buffer) Conc. 2; 10; 50; 200; 1000 and 5000 kU/l Preservative* <0.003%	Ready for use. Store at 2 – 8 °C until expiration date.

ImmunoCAP Total IgE Anti-IgE	
(mouse monoclonal antibodies) Preservative* <0.0015%	Ready for use. Store at 2 – 8 °C until expiration date.

Development Solution	
4-Methylumbelliferyl-β-D-galactoside 0.01% Preservative* <0.0010%	Ready for use. Store at 2 – 8 °C until expiration date. Do not freeze!

Stop Solution	
Sodium carbonate 4%	Ready for use. Store at 2 – 32 °C until expiration date.

Washing Solution	
For information, see separate Directions for Use for Washing Solution.	

ImmunoCAP Total IgE Control LMH	
Total IgE Control L (Low); Sodium azide 0.05%	Ready for use. Store at 2 – 8 °C until expiration date.
Total IgE Control M (Medium); Sodium azide 0.05%	
Total IgE Control H (High); Sodium azide 0.05%	
ImmunoCAP Total IgE Controls LMH are prepared from selected pooled human samples.	

ImmunoCAP IgE/ECP/Tryptase Sample Diluent	
(buffer solution with Bovine Serum Albumin) Preservative* <0.003%	Ready for use. Store at 2 – 8 °C until expiration date.

***Preservative:** Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1).

Additional material

- ImmunoCAP IgE/ECP/Tryptase Sample Diluent (10-9256-01/10-9360-01/10-9541-01)

Materials required but not provided by Phadia AB:

- Measuring cylinder 1000 ml
- Purified water (10, 11) or Clinical Laboratory Reagent Water (CLRW, 12)

⚠️ Precautions

- For in vitro diagnostic use. Not for internal or external use in humans or animals.
- Some reagents are manufactured from human blood components. The source materials have been tested by immunoassay for hepatitis B surface antigen, for antibodies to HIV1, HIV2 and hepatitis C virus and found to be negative. Nevertheless, all recommended precautions for the handling of blood derivatives should be observed. Please refer to Human Health Service (HHS) Publication No. (CDC) 93-8395 or other local/national guidelines on laboratory safety procedures.
- Reagents containing >0.0015% mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) may cause sensitization by skin contact. Avoid contact with skin. Wear suitable gloves. For more information see Safety Data Sheet.

INSTRUMENTS

Phadia 100 and Phadia 200 process all steps of the assay and print results automatically after the assay is completed^(a). Phadia 250, Phadia 1000, Phadia 2500 and Phadia 5000 are continuous random access systems that perform all steps of the assay^(a).

On board stability

Phadia 100 has no provisions for on board reagent storage. Information on reagent storage for the intended locations on the other Phadia instruments is shown below.

	Phadia 200	Phadia 250	Phadia 1000	Phadia 2500, Phadia 5000
Calibrator/ Curve Control		28 days.		N/A

	Phadia 200	Phadia 250	Phadia 1000	Phadia 2500, Phadia 5000
Conjugate	4 days. Recap bottles every night.		7 days.	
Development Solution	5 days. Recap bottles every night.		14 days.	5 days.
Stop Solution	14 days. Recap bottles every night.		14 days.	31 days.
Washing Solution (prepared solution)	7 days. Discard every seventh day and perform weekly maintenance according to respective instrument user manual.			
Washing Solution Concentrate	N/A			30 days.
Washing Solution Additive	N/A			30 days.
Sample Diluent	7 days. Recap bottles every night.		N/A	
ImmunoCAP Carrier	No provision for on board storage.	Until expiration date.		

SPECIMEN COLLECTION AND PREPARATION

Specimen collection
Serum and plasma (EDTA or heparin) samples from venous or capillary blood can be used. Collect blood samples using standard procedures. Keep specimens at room temperature (RT) for shipping purposes only. Store at 2 °C to 8 °C up to one week, or else at -20 °C. Avoid repeated freezing and thawing (9). For further reading on interfering substances, see reference (8).
Note: It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory. In general, laboratories should perform validation studies before implementing a change in specimen acceptance criteria (14).

Preparation of samples
Sample dilution is usually not required.
For determination of values higher than 5000 kU/l IgE, dilute the samples with ImmunoCAP IgE/ECP/Tryptase Sample Diluent.

Handling of control specimen
It is recommended to remove and recap the control vials from the instrument as soon as the pipetting of the samples is finished and the sample incubation is started. It is also recommended to gently stir the vial before use.
ImmunoCAP Total IgE Controls are ready for use and must not be further diluted. ImmunoCAP Total IgE Controls should be treated in the same way as a patient sample in the procedure.

PROCEDURES

Procedural steps^(a)
For procedural steps, see **Notes a**.

Parameters of the procedure
Patient samples are run in single determinations.
Volumes per determination:
Sample 40 µl
Conjugate 50 µl
Development Solution 50 µl
Stop Solution 600 µl

Incubations are performed at 37 °C by Phadia instruments.
Phadia 100: Total time for one assay is 2.5 hours.
Phadia 200: Total time for one assay is up to 4 hours.
Phadia 250, Phadia 1000, Phadia 2500 and Phadia 5000: The process time is 1 hour and 45 minutes from entering the first sample until the result is available.

Calibration^(a)
ImmunoCAP Total IgE Calibrators and ImmunoCAP Total IgE Calibrator Strip are run in duplicate to obtain a full calibration curve. The curve can be stored. The software for Phadia instruments has built-in acceptance limits for the calibration curve and the curve controls. Use two curve controls, CC-1 and CC-2, each in single determination to evaluate subsequent runs against the stored curve.
Calibrator range: 2-5000 kU/l.
Reference material: The IgE calibrators are traceable (via an unbroken chain of calibrations) to the 2nd International Reference Preparation (IRP) 75/502, or the equivalent 3rd International Standard 11/234, of Human Serum Immunoglobulin E from World Health Organization (WHO) (13).

QUALITY CONTROL

Record keeping for each assay
It is good laboratory practice to record the lot numbers of the components used, the dates when they were first opened and the remaining volumes.
Control specimen
Good laboratory practice requires that quality control specimen should be included in every run. Any material used should be assayed repeatedly to establish mean values and acceptable ranges.
Controls available from Phadia AB for day to day quality control:
• ImmunoCAP Total IgE Control LMH (10-9447-01)
Intended use
ImmunoCAP Total IgE Controls are used for monitoring ImmunoCAP Total IgE measurements performance in Phadia instruments.

Expected values for ImmunoCAP Total IgE Control LMH
As with all immunoassays the results are affected by the testing procedures and equipment used by different laboratories. It is therefore recommended that each laboratory establishes its own target value for each actual lot of control together with criteria of acceptance (recommended range ± 30%). This established target value is expected to fall within the range for the actual lot stated on the vial label. The range for each specific lot is calculated as a mean ± 3 SD using the expected long term variation. The mean value for each specific lot has been determined from 6 consecutive control assays, each in 6 replicates in ImmunoCAP Total IgE.

Proficiency testing
An external quality assessment program (proficiency testing) is available from Phadia AB for quality assurance purposes (Quality Club):
• Quality Club Total IgE (10-9297-01)

RESULTS

For Phadia instruments using Phadia Information Data Manager Software or Phadia Prime, all calculations are automatically performed^(a). Phadia 100 is programmed to automatically calculate all results^(a). It can also be used with Phadia Information Data Manager Software.

LIMITATIONS OF THE PROCEDURE

A definitive clinical diagnosis should only be made by the physician after all clinical and laboratory findings have been evaluated. It should not be based on the results of any single diagnostic method.

EXPECTED VALUES^(b)

Good laboratory practice recommends that each laboratory establishes its own expected range of values.
The results below have been obtained using Phadebas IgE PRIST. The expected values for Phadebas IgE PRIST can be used for ImmunoCAP Total IgE as comparison studies^(b) have shown good concordance between ImmunoCAP Total IgE and Pharmacia CAP System IgE FEIA, as well as between Pharmacia CAP System IgE FEIA and Phadebas PRIST.
Adults: Total IgE levels have been determined using Phadebas IgE PRIST in serum of 412 adult patients with respiratory symptoms, of which 160 were classified non-atopic and 252 had atopic disease, and showed the following distribution pattern between atopic and non-atopic individuals (7): below 25 kU/l - 84% non-atopic, above 100 kU/l - 78% atopic. When determined using Phadebas IgE PRIST (7), the geometric mean calculated from the total IgE levels in serum of 175 non-atopic adults was 13.2 kU/l, + 2 SD = 114 kU/l.
Children: The data from two independent studies using Phadebas IgE PRIST for the determination of total IgE in serum of 466 carefully selected healthy children (4,5), have been used for calculations leading to the following summary of development of serum total IgE levels during childhood (6). After the peak at the age of 10 years, serum total IgE levels decline to adult values.

Age	Geometric mean (kU IgE/l)	+1 SD (kU IgE/l)
6 weeks	0.6	2.3
3 months	1.0	4.1
6 months	1.8	7.3
9 months	2.6	10
12 months	3.2	13
2 years	5.7	23
3 years	8.0	32
4 years	10	40
5 years	12	48
6 years	14	56
7 years	16	63
8 years	18	71
9 years	20	78
10 years	22	85

PERFORMANCE CHARACTERISTICS

Instrument comparison
Comparison studies^(b) have been performed with different combinations of Phadia instruments, including more than 85 patient samples. The results obtained show good concordance between instruments.

Precision^(b)
For Phadia 100
The following pooled coefficients of variation have been obtained. Each sample has been assayed in 4 replicates on 18 different occasions, using stored calibration curves.

Sample level (kU/l)	Coefficients of variation (%)	
	Within assay	Between assay
13	3	9
75	2	5
640	3	9

For Phadia 200, Phadia 250, Phadia 1000, Phadia 2500 and Phadia 5000
The following pooled coefficients of variation have been obtained with Phadia 1000. Each sample has been assayed in 2 replicates on 24 different occasions, using the same lot of reagents. Each level contains 3 samples. The values are also representative for Phadia 200, Phadia 250, Phadia 2500 and Phadia 5000.

Sample level (kU/l)	Coefficients of variation (%)	
	Within assay	Between assay
15-60	3	5
75-430	3	4
600-1840	3	7

Analytical sensitivity^(b)
The detection limit is 2 kU/l.












Analytical specificity^(b)
The cross-reactivity with other human immunoglobulins is non-detectable at physiological concentrations of IgA, IgD, IgM and IgG.

Recovery^(b)
Mean recovery is 98%.

WARRANTY

The performance data presented here was obtained using the procedure indicated. Any change or modification in the procedure not recommended by Phadia AB may affect the results, in which event Phadia AB disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. Phadia AB and its authorized distributors, in such event, shall not be liable for damages indirect or consequential.

SYMBOLS

	Use-by date		Contains sufficient for <n> tests
	Batch code		In vitro diagnostic medical device
	Date of manufacture		Temperature limit
	Catalogue number		Consult instructions for use
	Caution		Biological risks
	Manufacturer		

Full symbol glossary is available at: https://symbols_glossary.phadia.com.

REFERENCES

- Ishizaka K, Ishizaka T, Hornbrook MM. J Immunol 1966;97(1);75-85.
- Johansson SGO, Bennich HH, Wide L. A new class of immunoglobulin in human serum. Immunology 1968;14:265-272.
- Bennich H, Ishizaka K, Johansson SGO, Rowe DS, Stanworth DR, Terry WD. WHO Bull 1968;38;151-2.
- Bhalla RB, Rappaport I, De Filippi I, Schwartz MK. Serum IgE Levels in a northeast United States caucasian population. In: Heusghem C, Albert A, Benson E S, eds. Advanced interpretation of clinical laboratory data. New York:Marcel Dekker Inc, 1982:295-305.
- Kjellman NI, Johansson SGO, Roth A. Clin Allergy 1976;6:51-9.
- Björkstén B, Weeke B. Allergy 1985;40 (Suppl 4).
- Zetterström O, Johansson SGO. IgE concentrations measured by PRIST in serum of healthy adults and in patients with respiratory allergy. Allergy 1981;36:537-47.
- Tate J, Ward G. Interference in Immunoassay. Clin Biochem Rev 2004;25:105-120.
- CSLI I/LA20-A2. Analytical Performance Characteristics and Clinical Utility of Immunological Assays for Human Immunoglobulin E (IgE) Antibodies and Defined Allergen Specificities; Approved Guideline – Second Edition, 2009.
- US Pharmacopeia & National Formulary, current edition.
- European Pharmacopeia, current edition.
- CLSI C3-A4. Preparation and Testing of Reagent Water in the Clinical Laboratory, Approved Guideline-Fourth Edition, Clinical and Laboratory Standards Institute, 2006.
- Thorpe S, Heath A, Fox B, Patel D, Egner W. The 3rd International Standard for serum IgE: international collaborative study to evaluate a candidate preparation. Clin Chem Lab Med 2014; 52(9); 1283-1289.
- CLSI. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline-Forth Edition. CLSI document H18-A4, Wayne, PA: Clinical and Laboratory Standards Institute; 2010.

Notes

^(a)For more information, see Phadia 100, Phadia 200, Phadia 250, Phadia 1000 and/or Phadia 2500/Phadia 5000 User Manual.

^(b)Studies performed at Phadia AB, Uppsala, Sweden.

Patents/Trademarks

The following designations are trademarks belonging to Phadia AB: ImmunoCAP, Phadia, Phadia, Quality Club.

Trademark change: Phadia AB has changed the trademarks of the instrument platforms from "UniCAP™" and "ImmunoCAP™" to "Phadia™". The new name has been applied to the instruments and related items, e.g. Software and User Manuals. The trademark "ImmunoCAP™" has been removed from the System Reagents. This is a trademark change only; the change has no impact on performance or safety.

Addresses

AUSTRIA Thermo Fisher Diagnostics Austria GmbH
Dresdner Str. 89
A-1200 Vienna
Tel: +43-1 270 2020 Fax: +43-1 270 202020
BELGIUM Thermo Fisher Diagnostics NV
Pontbeekstraat 2
BE-1702 GROOT-BIJGAARDEN
Tel: +32-2 749 55 15 Fax: +32-2 749 55 23
BRAZIL Phadia Diagnósticos Ltda.
Rua Eugênio de Medeiros, 303 cj 1101C
05425-000 São Paulo – SP
Tel: +55 0800 5515 355
CHINA Thermo Fisher Scientific (China) Co., Ltd.
Building 6-7
No. 27 Xin Jin Qiao Road
Shanghai 201102
P.R. China
Tel: +86 800 810 5118 Fax: +86 21 6100 21288
CZECH REPUBLIC
Thermo Fisher Scientific, Phadia s.r.o.
Drahobejlova 1019/27
19000 PRAHA 9
Tel: +420 220 518 743 Fax: +420 220 518 743
DENMARK Thermo Fisher Diagnostics ApS
Gydevang 33
DK-3450 ALLERØD
Tel: +45-70 23 33 06 Fax: +45-70 23 33 07
FINLAND Phadia Oy
Ratastie 2
P.O. Box 100
FI-01621
Tel: +358 10 3292 110 Fax: +358 10 3292 531
FRANCE Thermo Fisher Diagnostics S.A.S.
CS60610
FR-78056 ST QUENTIN-YVELINES CEDEX
Tel: +33-1 61 37 34 30 Fax: +33-1 30 64 62 37
GERMANY Thermo Fisher Diagnostics GmbH
Munzinger Str. 7, D-79111 FREIBURG
Tel: +49-761 47 805-0 Fax: +49-761 47805-338
Mail: Thermo Fisher Diagnostics GmbH
Postfach 1050 DE-790 10 FREIBURG
HONG KONG
Thermo Fisher Scientific (Hong Kong) Limited
Units 01, 09-13 & 15-18, Level 13, Tower 1
Kowloon Commerce Centre
No. 51 Kwai Cheong Road, Kwai Chung
New Territories, Hong Kong
Tel: +852 3107 7600 Fax: +852 2567 4447
INDIA Thermo Fisher Scientific India Pvt. Ltd
Unit No.07, 10 & 11, ground floor,
Splendor forum, plot no 03
Jasola, District Centre,
NEW DEHLI-110025
Tel: +91 11 49375400
IRELAND Thermo Fisher Diagnostics (Ireland) Ltd.
16 Shenley Pavilions, Chalkdell Drive
Shenley Wood, Milton Keynes, MK5 6LB
Tel: +44 1908 76 91 10 Fax: +44 844 324 94
ITALY Thermo Fisher Diagnostics S.p.A.
Via G. B. Tiepolo, 18
IT-20900 MONZA (MB)
Tel: +39 039 8389.1 Fax: +39 039 838 9329
JAPAN Thermo Fisher Diagnostics K.K.
Sumitomo Fudosan Mita Twin Bldg. East 4F
4-2-8 Shibaura, Minato-ku,
Tokyo 108-0023, Japan
Tel: +81 3 6872 6200 Fax: +81 3 6872 6220
KOREA Thermo Fisher Scientific Korea Ltd.
11-12FL, Suseo Office Building
281, Gwangpyeong-ro, Gangnam-gu,
Seoul 06349, Korea
Tel: +82 2 6196 5556-9 Fax: +82 2 6196 5555
THE NETHERLANDS Thermo Fisher Diagnostics B.V.
Postbus 696
NL-3430 AR NIEUWEGEIN

Tel: +31-30 602 37 00 Fax: +31-30 602 37 09

NORWAY Thermo Fisher Diagnostics AS
Postboks 114, Smedstad
NO-0309 OSLO

Tel: +47-21 67 32 80

PORTUGAL Thermo Fisher Diagnostics,
Sociedade Unipessoal Lda
Lagoas Park - Edifício nº11 - Piso 0
PT-2740-270 PORTO SALVO

Tel: +351-214 23 53 50 Fax: +351-214 21 60 36

SOUTH AFRICA Laboratory Specialities (Pty)
t/a Thermo Fisher Scientific

Thermo Fisher Scientific Building

Unit A, 3 Susan Street

Strijdom Park, Randburg 2169

Tel: +27 11 792 6790 Fax: +27 11 793 1064

SPAIN Thermo Fisher Diagnostics S.L.U.

Avda. Alcalde Barnils nº 70 Planta 2

Edificio Onada

08174 Sant Cugat del Vallés, Barcelona

Tel: +34-935 765 800 Fax: +34-935 765 820

SWEDEN Thermo Fisher Diagnostics AB

C/o Phadia AB

P O Box 6460 SE-751 37 UPPSALA

Tel: +46-18 16 60 60 Fax: +46-18 16 63 24

SWITZERLAND Thermo Fisher Diagnostics AG

Sennweidstrasse 46

CH-6312 STEINHAUSEN

Tel: +41-43 343 40 50 Fax: +41-43 343 40 51

TAIWAN Phadia Taiwan Inc.

6F-1, No. 85, Zhouzih St., NeiHu District

Taipei, Taiwan, Republic of China

Tel: +886 2 8751 6655 Fax: +886 2 8751 5353

UNITED KINGDOM Thermo Fisher Diagnostics Ltd

16 Shenley Pavilions, Chalkdell Drive

Shenley Wood, Milton Keynes, MK5 6LB

Tel: +44-1908 76 91 10 Fax: +44-844 324 94 95

USA Phadia US Inc.

4169 Commercial Avenue

Portage, Michigan 49002

Tel: +1 800-346-4364 (Toll Free) Fax: +1 269 492-7541

OTHER COUNTRIES Phadia AB,

Distributor Sales

P O Box 6460, SE-751 37 UPPSALA

Tel: +46 18 16 50 00 Fax: +46 18 16 63 65



Phadia AB,
Rapskatan 7P, P. O. Box 6460, 751 37 Uppsala, Sweden
Tel: +46 18 16 50 00 Fax: +46 18 14 03 58



Issued February 2012. Revised September 2018.

© Phadia AB, Uppsala, Sweden.