

ImmunoCAP™ Tryptase

Fluoroenzymeimmunoassay

Directions for Use 52-5467-EN/03

INTENDED USE

ImmunoCAP Tryptase is an in vitro test system for the quantitative measurement of tryptase in human serum or plasma. It is intended for in vitro diagnostic use in conjunction with other clinical findings, and is to be used in clinical laboratories. ImmunoCAP Tryptase is to be used with the instruments Phadia 100, Phadia 200, Phadia 250 or Phadia 1000.

SUMMARY AND EXPLANATION OF THE TEST

Tryptase is the most abundant protein in mast cells. During IgE mediated allergic reactions mast cells are activated and release inflammatory mediators including tryptase (1,2,3). ImmunoCAP Tryptase measures the total tryptase levels including all forms of α -tryptase and β -tryptase.

The baseline level of tryptase in the circulation reflects the number of mast cells.

Elevated basal levels of serum tryptase and/or an underlying mastocytosis may be risk factors particularly in patients with history of severe reactions. This should for example be taken into consideration in venom immunotherapy (8-10).

Mature β -tryptase is transiently elevated in most cases of systemic anaphylactic reactions. The peak level is usually reached 15-120 minutes after onset of the reaction and tryptase levels then decline slowly within the next 3-6 hours. The return to baseline levels can generally be verified approximately 24 hours after the reaction (1,7).

Elevated tryptase levels in post-mortem samples may indicate a fatal anaphylactic reaction as a cause of death (11,12).

PRINCIPLE OF THE PROCEDURE

Anti-tryptase, covalently coupled to ImmunoCAP, reacts with the tryptase in the patient sample. After washing, enzyme labeled antibodies against tryptase are added to form a complex. Following incubation, unbound enzyme-anti-tryptase is washed away and the bound complex is then incubated with a developing agent. After stopping the reaction, the fluorescence of the eluate is measured. The higher the response value, the more tryptase is present in the sample. To evaluate the test results, the responses for the patient samples are transformed to concentrations with the use of a calibration curve.

REAGENTS AND MATERIAL

Reagents are packaged as described below, each purchased separately. The two digit suffix (-xx) on the article number may vary between countries. All kits are not available in all countries. The expiration date and storage temperature are stated on the labels. Do not use reagents beyond their expiration dates.

Note: It is not recommended to pool any reagents.

Keep the ImmunoCAP carrier closed to avoid evaporation of buffer.

Reagents for Phadia 100

- **ImmunoCAP Tryptase** (Art No 10-9525-01: for 48 determinations)
 - Tryptase Conjugate (1 vial)
 - Tryptase Curve Control 1 (CC-1) (4 single dose vials)
- **ImmunoCAP Tryptase Calibrators** (Cal-xx) (Art No 10-9526-01: for 1 calibration curve)
- **ImmunoCAP Tryptase Curve Control** (CC-1) (Art No 10-9527-01: for 6 additional assay runs)
- **ImmunoCAP Tryptase Anti-Tryptase** (a-Tryp) (Art No 14-4518-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9478-01: for 600 determinations)
- **Stop Solution** (Art No 10-9479-01: for 600 determinations)
- **Washing Solution** (Art No 10-9422-01: 6 x 1 l)
 - Washing Solution Additive, 6 x 17.2 ml
 - Washing Solution Concentrate, 6 x 80 ml
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml
- **ImmunoCAP Tryptase Control** (Art No 10-9370-01: 6 vials)

Reagents for Phadia 200

- **ImmunoCAP Tryptase Conjugate 50** (Art No 10-9522-01: for 2 x 50 determinations)
- **ImmunoCAP Tryptase Calibrator Strip** (Art No 10-9523-01: for 1 calibration curve)
- **ImmunoCAP Tryptase Curve Control Strip** (CC-1) (Art No 10-9524-01: for 6 x 1 curve controls)
- **ImmunoCAP Tryptase Anti-Tryptase** (a-Tryp) (Art No 14-4518-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9441-01: for 6 x 200 determinations; Art No 10-9440-01: for 6 x 315 determinations)

- **Stop Solution** (Art No 10-9442-01: for 6 x 185 determinations; Art No 10-9479-01: for 6 x 100 determinations)
- **Washing Solution** (Art No 10-9422-01: 6 x 1 l)
 - Washing Solution Additive, 6 x 17.2 ml
 - Washing Solution Concentrate, 6 x 80 ml
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml
- **ImmunoCAP Tryptase Control** (Art No 10-9370-01: 6 vials)

Reagents for Phadia 250

- **ImmunoCAP Tryptase Conjugate 50** (Art No 10-9522-01: for 2 x 50 determinations)
- **ImmunoCAP Tryptase Calibrator Strip** (Art No 10-9523-01: for 1 calibration curve)
- **ImmunoCAP Tryptase Curve Control Strip** (CC-1) (Art No 10-9524-01: for 6 x 1 curve controls)
- **ImmunoCAP Tryptase Anti-Tryptase** (a-Tryp) (Art No 14-4518-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9441-01: for 6 x 200 determinations; Art No 10-9440-01: for 6 x 315 determinations)
- **Stop Solution** (Art No 10-9442-01: for 6 x 185 determinations)
- **Washing Solution** (Art No 10-9422-01: 6 x 1 l)
 - Washing Solution Additive, 6 x 17.2 ml
 - Washing Solution Concentrate, 6 x 80 ml
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml
- **ImmunoCAP Tryptase Control** (Art No 10-9370-01: 6 vials)

Reagents for Phadia 1000

- **ImmunoCAP Tryptase Conjugate 50** (Art No 10-9522-01: for 2 x 50 determinations)
- **ImmunoCAP Tryptase Calibrator Strip** (Art No 10-9523-01: for 1 calibration curve)
- **ImmunoCAP Tryptase Curve Control Strip** (CC-1) (Art No 10-9524-01: for 6 x 1 curve controls)
- **ImmunoCAP Tryptase Anti-Tryptase** (a-Tryp) (Art No 14-4518-01: carrier of 16 ImmunoCAP)
- **Development Solution** (Art No 10-9439-01: for 6 x 1200 determinations; Art No 10-9314-01: for 6 x 2000 determinations)
- **Stop Solution** (Art No 34-2271-51: for 1200 determinations)
- **Washing Solution** (Art No 10-9202-01: 2 x 5 l)
 - Washing Solution Additive, 2 x 86 ml
 - Washing Solution Concentrate, 2 x 400 ml
- **ImmunoCAP Tryptase Control** (Art No 10-9370-01: 6 vials)

Details of reagents

ImmunoCAP Tryptase Conjugate	
β -Galactosidase-anti-tryptase Approximately 1.6 μ g/ml (mouse monoclonal antibodies) Sodium azide 0.06%	Ready for use. Store at 2 – 8 °C until expiration date. Do not freeze!

ImmunoCAP Tryptase Calibrators/Calibrator Strip	
(human recombinant tryptase in buffer) Conc. 1; 5; 12.5; 50 and 200 μ g/l Sodium azide 0.05%	Ready for use. Store at 2 – 8 °C until expiration date.

ImmunoCAP Tryptase Curve Control/Curve Control Strip	
(human recombinant tryptase in buffer) Sodium azide 0.05%	Ready for use. Store at 2 – 8 °C until expiration date.

ImmunoCAP Tryptase Anti-Tryptase	
(mouse monoclonal antibodies) Preservative* <0.0015%	Ready for use. Store at 2 – 8 °C until expiration date.

Development Solution	
4-Methylumbelliferyl- β -D-galactoside 0.01% Preservative* <0.0010%	Ready for use. Store at 2 – 8 °C until expiration date. Do not freeze!

Stop Solution	
Sodium carbonate 4%	Ready for use. Store at 2 – 32 °C until expiration date.

Washing Solution	
For information, see separate Directions for Use for Washing Solution.	

ImmunoCAP Tryptase Control	
ImmunoCAP Tryptase Control is prepared from selected pooled human sera. It is lyophilized to ensure maximum stability. Reconstitution of ImmunoCAP Tryptase Control Reconstitute the content of a vial by adding exactly 500 μ l purified water (18, 19) or clinical laboratory reagent water (CLRW, 20). Let the vial stand for one minute, then mix gently until the content is completely dissolved. Shelf-life and storage Lyophilized serum: Store at 2-8 °C until expiration date. Reconstituted serum: Store at 2 – 8 °C for 1 week or -20 °C for 4 weeks. (Repeated freezing and thawing should be avoided.)	

ImmunoCAP IgE/ECP/Tryptase Sample Diluent	
(buffer solution with Bovine Serum Albumin) Preservative* <0.003%	Ready for use. Store at 2 – 8 °C until expiration date.

***Preservative:** Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1).

Additional material

Additional products available from Phadia AB:

- ImmunoCAP IgE/ECP/Tryptase Sample Diluent (10-9256-01/10-9360-01/10-9541-01)

Materials required but not provided by Phadia AB:

- Measuring cylinder 1000 ml
- Purified water (18,19) or Clinical Laboratory Reagent Water (CLRW, 20)



Precautions

- For in vitro diagnostic use. Not for internal or external use in humans or animals.
- Some reagents are manufactured from human blood components. The source materials have been tested by immunoassay for hepatitis B surface antigen, for antibodies to HIV1, HIV2 and hepatitis C virus and found to be negative. Nevertheless, all recommended precautions for the handling of blood derivatives should be observed. Please refer to Human Health Service (HHS) Publication No. (CDC) 93-8395 or other local/national guidelines on laboratory safety procedures.
- Reagents containing >0.0015% mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) may cause sensitization by skin contact. Avoid contact with skin. Wear suitable gloves. For more information see Safety Data Sheet.

INSTRUMENTS

Phadia 100 and Phadia 200 process all steps of the assay and print results automatically after the assay is completed^(a). Phadia 250 and Phadia 1000 are continuous random access systems that perform all steps of the assay^(a).

On board stability

Phadia 100 has no provisions for on board reagent storage. Information on reagent storage for the intended locations on the other Phadia instruments is shown below.

	Phadia 200	Phadia 250	Phadia 1000
Calibrator/ Curve Control	28 days.		
Conjugate	4 days. Recap bottles every night.		7 days.
Development Solution	5 days. Recap bottles every night.		14 days.

	Phadia 200	Phadia 250	Phadia 1000
Stop Solution	14 days. Recap bottles every night.		14 days.
Washing Solution (prepared solution)	7 days. Discard every seventh day and perform weekly maintenance according to respective instrument user manual.		
Washing Solution Concentrate	N/A		
Washing Solution Additive	N/A		
Sample Diluent	7 days. Recap bottles every night.		N/A
ImmunoCAP Carrier	No provision for on board storage.	Until expiration date.	

SPECIMEN COLLECTION AND PREPARATION

Specimen collection

Serum and plasma (EDTA or heparin) samples from venous blood can be used. Collect blood samples and prepare serum or plasma according to standard procedures. Keep specimens at room temperature (RT) for shipping purposes only, up to 2 days. Store at 2 °C to 8 °C for up to one week, or else at –20 °C. Avoid repeated freezing and thawing (16).

It is preferable that samples be taken no earlier than 15 minutes from onset / up to 3 hours after the onset of the suspected incident causing mast cell activation (5,6). The time between the reaction and sample collection should be noted. To confirm the return to baseline levels an additional blood sample should be collected after 24 - 48 hours, time depending on the magnitude of the activation. At the suspicion of elevated basal levels or an underlying mastocytosis additional sample(s) should be taken 1-2 week(s) later. Post mortem samples should be taken within 48 hours from time of death.

Note: It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory. In general, laboratories should perform validation studies before implementing a change in specimen acceptance criteria (21).

Preparation of samples

Sample dilution is usually not required.

For determination of values higher than 200 µg/l, dilute the samples with ImmunoCAP IgE/ECP/Tryptase Sample Diluent.

Handling of control specimen

It is recommended to remove and recap the control vials from the instrument as soon as the pipetting of the samples is finished and the sample incubation is started. It is also recommended to gently stir the vial before use. Avoid repeated freezing and thawing. ImmunoCAP Tryptase Control should be treated in the same way as a patient sample in the procedure.

PROCEDURES

Procedural steps^(a)

For procedural steps, see **Notes a**.

Parameters of the procedure

Patient samples are run in single determinations.

Volumes per determination:

Sample	40 µl
Conjugate	50 µl
Development Solution	50 µl
Stop Solution	600 µl

Incubations are performed at 37 °C by Phadia instruments.

Phadia 100: Total time for one assay is 2.5 hours.

Phadia 200: Total time for one assay is up to 4 hours.

Phadia 250 and Phadia 1000: The process time is 1 hour and 45 minutes from entering the first sample until the result is available.

Calibration^(a)

The ImmunoCAP Tryptase Calibrators or ImmunoCAP Tryptase Calibrator Strip are run in duplicates to obtain a full calibration curve. The curve can be stored. The software for Phadia instruments have built-in acceptance limits for the calibration curve and the curve controls.

Use one curve control, CC-1, in duplicate determinations to evaluate subsequent assays against the stored curve.

Calibrator range: 1-200 µg/l.

Reference material: ImmunoCAP Tryptase Calibrators are calibrated against an in-house tryptase reference which traces back to tryptase purified according to Schwartz et al (13).

QUALITY CONTROL

Record keeping for each assay

It is good laboratory practice to record the lot numbers of the components used, the dates when they were first opened and the remaining volumes.

Control specimen

Good laboratory practice requires that quality control specimen should be included in every run. Any material used should be assayed repeatedly to establish mean values and acceptable ranges.

Controls available from Phadia AB for day to day quality control:

- ImmunoCAP Tryptase Control (10-9370-01)

Intended use

ImmunoCAP Tryptase Control is used for monitoring ImmunoCAP Tryptase measurements performance in Phadia instruments.

Expected values for ImmunoCAP Tryptase Control

As with all immunoassays the results are affected by the testing procedures and equipment used by different laboratories. It is therefore recommended that each laboratory establishes its own target value for each actual lot of control together with criteria of acceptance (recommended range \pm 30 %).

This established target value is expected to fall within the range for the actual lot stated on the vial label.

The range for each specific lot is calculated as a mean \pm 3 SD using the expected long term variation. The mean value for each specific lot has been determined from 6 consecutive control assays, each in 6 replicates using ImmunoCAP Tryptase.

RESULTS

For Phadia instruments using Phadia Information Data Manager Software or Phadia Prime, all calculations are automatically performed^(a). Phadia 100 is programmed to automatically calculate all results^(a). It can also be used with Phadia Information Data Manager Software.

LIMITATIONS OF THE PROCEDURE

A definitive clinical diagnosis should only be made by the physician after all clinical and laboratory findings have been evaluated. It should not be based on the results of any single diagnostic method.

- Systemic anaphylactic reactions without elevation of tryptase might be caused by other non-mast cell pathways.
- In general, elevated tryptase levels are more commonly seen after parenteral than oral introduction of agents causing systemic anaphylactic reactions.
- Moderately elevated tryptase levels have been observed in post-mortem samples by other death causes than anaphylactic reactions.
- Heterophilic antibodies, especially human-anti mouse antibodies (HAMA), in human serum/plasma might react with the mouse immunoglobulins used as capturing antibodies on the solid phase and as detection antibodies in the conjugate in ImmunoCAP Tryptase (14,15). The presence of heterophilic antibodies is uncommon but can cause false results, mostly falsely elevated levels. The composition of the anti-tryptase conjugate is designed in order to minimize this kind of interference.

Still uncommon, the risk for interference is increased in certain patient groups e.g. patients having rheumatoid factor (RF) or patients receiving preparations containing mouse monoclonal antibodies (including chimeric/humanized) for diagnostic and/or therapeutic use. Another risk group is patients regularly exposed to animals and/or animal products.

One method to confirm the presence of heterophilic antibodies is by pre-treatment of samples using commercially available Heterophilic Blocking Tubes (HBT). For use and limitations of HBT, see manufacturer's instructions.

EXPECTED VALUES^(b)

Good laboratory practice recommends that each laboratory establishes its own expected range of values.

A study^(b) with 124 self-reported healthy individuals (56 males and 68 females) was performed with Phadia 250. The age range was 3-67 years for males and 4-63 years for females.

The following ImmunoCAP Tryptase results were obtained:

Geometric mean	3.4 µg/l
95 upper percentile	11.0 µg/l

In patients with systemic mastocytosis levels of tryptase are, in general, persistently elevated above 20 µg/l (4).

Baseline tryptase levels in the range of approximately 10-20 µg/l reflect an increased mast cell burden indicating an increased risk in patients with history of severe anaphylactic reaction (1).

In severe cases the triggering agent causing a transiently elevation of tryptase should be identified.

PERFORMANCE CHARACTERISTICS

Instrument comparison

Comparison studies^(b) have been performed with different combinations of Phadia instruments, including at least 50 patient samples. The results obtained show good concordance between instruments.

Precision^(b)

The following pooled coefficients of variation have been obtained with Phadia 250 and 4-11 samples/level. Each sample has been assayed in 4 replicates on 7 different occasions on each of 3 different instruments, using the same lot of reagents and stored calibration curves. The values are also representative for Phadia 100, Phadia 200 and Phadia 1000.

Sample level (µg/l)	Coefficients of variation (%)	
	Within assay	Between assay
1-20	3	5
20-100	3	6
100-200	4	7

Analytical sensitivity^(b)

The limit of detection (17) is 1.0 µg/l.

Analytical specificity^(b)

No interference of Heparin, Rheumatoid factor, hemolysed, lipemic or icteric samples has been observed.












Recovery^(b)

Mean recovery is 100 \pm 7.5%.

WARRANTY

The performance data presented here was obtained using the procedure indicated. Any change or modification in the procedure not recommended by Phadia AB may affect the results, in which event Phadia AB disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. Phadia AB and its authorized distributors, in such event, shall not be liable for damages indirect or consequential.

SYMBOLS

	Use-by date		Contains sufficient for <n> tests
	Batch code		In vitro diagnostic medical device
	Date of manufacture		Temperature limit
	Catalogue number		Consult instructions for use
	Caution		Biological risks
	Manufacturer		

Full symbol glossary is available at: https://symbols_glossary.phadia.com.

REFERENCES

1. Schwartz LB: Diagnostic Value of Tryptase in Anaphylaxis and Mastocytosis. *Immunol Allergy Clin N Am* 2006;26:451–463.
2. Schwartz LB: Mast Cells and Basophils. *Clin Allergy Immunol* 2002;16:3-42.
3. Caughey GH: Tryptase genetics and anaphylaxis. *J Allergy Clin Immunol* 2006;117(6):1411-1414.

4. Valent P: Diagnostic Evaluation and Classification of Mastocytosis. *Immunol Allergy Clin N Am* 2006;26:515-534.
5. Lieberman P, Nicklas RA, Oppenheimer J et al.: The Diagnosis and Management of Anaphylaxis Practice Parameter: 2010 Update. *J Allergy Clin Immunol* 2010 Sep;126(3):477-80.e1-42.
6. Ebo DG, Fisher MM, Hagendorens MM, Bridts CH, Stevens WJ: Anaphylaxis during anaesthesia: diagnostic approach. *Allergy* 2007;62:471-487.
7. Schwartz LB, Yunginger JW, Miller J, Bokhari R and Dull D: Time Course of Appearance and Disappearance of Human Mast Cell Tryptase in the Circulation after Anaphylaxis. *J Clin Invest* 1989;83:1551-1555.
8. Haeberli G, Brönnimann M, Hunziker T and Müller U: Elevated basal serum tryptase and hymenoptera venom allergy: relation to severity of sting reactions and to safety and efficacy of venom immunotherapy. *Clin Exp Allergy* 2003;33:1216-1220.
9. Biló BM, Rueff F, Mosbech H, Bonifazi F, Oude-Elberink JNG & the EAACI Interest Group on Insect Venom Hypersensitivity: Diagnosis of Hymenoptera venom allergy. *Allergy* 2005;60:1339-1349 / EAACI Position Paper <http://www.eaaci.net/media/PDF/D/652.pdf>.
10. Bonifazi F, Jutel M, Biló BM, Birnbaum J, Müller U and the EAACI Interest Group on Insect Venom Hypersensitivity: Prevention and treatment of hymenoptera venom allergy: guidelines for clinical practice. *Allergy* 2005;60:1459-1470/EAACI Position Paper <http://www.eaaci.net/media/PDF/P/653.pdf>.
11. Yunginger JW, Nelson DR, Squillace DL, Jones RT, Holley KE, Hyma BA, Biedrzycki L, Sweeney KG, Stumer WQ and Schwartz LB: Laboratory Investigation of Deaths Due to Anaphylaxis. *J Forensic Sci* 1991;36:857-865.
12. Edston E, van Hage-Hamsten M: Mast cell tryptase and hemolysis after trauma. *Forensic Science International* 2003;131:8-13.
13. Schwartz L.B., et al.: Immunologic and physicochemical evidence for conformational changes occurring on conversion of human mast cell tryptase from active tetramer to inactive monomer. *J Immunol.* 1990 (144) 2304-2311.
14. Bjerner J, Børner OP, Nustad K: The War on Heterophilic Antibody Interference. *Clin Chem* 2005;51(1):9-11.
15. Kricka LJ: Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999;45(7):942-956.
16. CLSI H18-A4. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline - Fourth Edition, Clinical and Laboratory Standards Institute, 2010.
17. NCCLS EP17-A. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline – First Edition, National Committee for Clinical Laboratory Standards, 2004.
18. US Pharmacopeia & National Formulary, current edition.
19. European Pharmacopeia, current edition.
20. CLSI C3-A4, Preparation and Testing of Reagent Water in the Clinical Laboratory, Approved Guideline-Fourth Edition, Clinical and Laboratory Standards Institute, 2006.
21. CLSI. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline-Forth Edition. CLSI document H18-A4, Wayne, PA: Clinical and Laboratory Standards Institute; 2010.

Notes

(a)For more information, see Phadia 100, Phadia 200, Phadia 250 and/or Phadia 1000 User Manual.

(b)Studies performed at Phadia AB, Uppsala, Sweden.

Patents/Trademarks

The following designations are trademarks belonging to Phadia AB: ImmunoCAP, Phadia.

Trademark change: Phadia AB has changed the trademarks of the instrument platforms from "UniCAP™" and "ImmunoCAP™" to "Phadia™". The new name has been applied to the instruments and related items, e.g. Software and User Manuals. The trademark "ImmunoCAP™" has been removed from the System Reagents. This is a trademark change only; the change has no impact on performance or safety.

Addresses

AUSTRIA Thermo Fisher Diagnostics Austria GmbH
Dresdner Str. 89
A-1200 Vienna
Tel: +43-1 270 2020 Fax: +43-1 270 202020
BELGIUM Thermo Fisher Diagnostics NV
Pontbeekstraat 2
BE-1702 GROOT-BIJGAARDEN
Tel: +32-2 749 55 15 Fax: +32-2 749 55 23
BRAZIL Phadia Diagnósticos Ltda.
Rua Eugênio de Medeiros, 303 cj 1101C
05425-000 São Paulo – SP
Tel: +55 0800 5515 355

CHINA Thermo Fisher Scientific (China) Co., Ltd.
Building 6-7

No. 27 Xin Jin Qiao Road
Shanghai 201102
P.R. China

Tel: +86 800 810 5118 Fax: +86 21 6100 21288

CZECH REPUBLIC

Thermo Fisher Scientific, Phadia s.r.o.

Drahobějova 1019/27

19000 PRAHA 9

Tel: +420 220 518 743 Fax: +420 220 518 743

DENMARK

Thermo Fisher Diagnostics ApS

Gydevang 33

DK-3450 ALLERØD

Tel: +45-70 23 33 06 Fax: +45-70 23 33 07

FINLAND

Phadia Oy

Ratatie 2

P.O. Box 100

FI-01621

Tel: +358 10 3292 110 Fax: +358 10 3292 531

FRANCE

Thermo Fisher Diagnostics S.A.S.

CS60610

FR-78056 ST QUENTIN-YVELINES CEDEX

Tel: +33-1 61 37 34 30 Fax: +33-1 30 64 62 37

GERMANY

Thermo Fisher Diagnostics GmbH

Munzinger Str. 7, D-79111 FREIBURG

Tel: +49-761 47 805-0 Fax: +49-761 47805-338

Mail: Thermo Fisher Diagnostics GmbH

Postfach 1050 DE-790 10 FREIBURG

HONG KONG

Thermo Fisher Scientific (Hong Kong) Limited

Units 01, 09-13 & 15-18, Level 13, Tower 1

Kowloon Commerce Centre

No. 51 Kwai Cheong Road, Kwai Chung

New Territories, Hong Kong

Tel: +852 3107 7600 Fax: +852 2567 4447

INDIA

Thermo Fisher Scientific India Pvt. Ltd

Unit No.07, 10 & 11, ground floor,

Splendor forum, plot no 03

Jasola, District Centre,

NEW DEHLI-110025

Tel: +91 11 49375400

IRELAND

Thermo Fisher Diagnostics (Ireland) Ltd.

16 Shenley Pavilions, Chalkdell Drive

Shenley Wood, Milton Keynes, MK5 6LB

Tel: +44 1908 76 91 10 Fax: +44 844 324 94

ITALY

Thermo Fisher Diagnostics S.p.A.

Via G. B. Tiepolo, 18

IT-20900 MONZA (MB)

Tel: +39 039 8389.1 Fax: +39 039 838 9329

JAPAN

Thermo Fisher Diagnostics K.K.

Sumitomo Fudosan Mita Twin Bldg. East 4F

4-2-8 Shibaura, Minato-ku,

Tokyo 108-0023, Japan

Tel: +81 3 6872 6200 Fax: +81 3 6872 6220

KOREA

Thermo Fisher Scientific Korea Ltd.

11-12FL, Suseo Office Building

281, Gwanggyeong-ro, Gangnam-gu,

Seoul 06349, Korea

Tel: +82 2 6196 5556-9 Fax: +82 2 6196 5555

THE NETHERLANDS

Thermo Fisher Diagnostics B.V.

Postbus 696

NL-3430 AR NIEUWEGEIN

Tel: +31-30 602 37 00 Fax: +31-30 602 37 00

NORWAY

Thermo Fisher Diagnostics AS

Postboks 114, Smestad

NO-0309 OSLO

Tel: +47-21 67 32 80

PORTUGAL

Thermo Fisher Diagnostics,

Sociedade Unipessoal Lda

Lagoas Park - Edifício n°11 - Piso 0

PT-2740-270 PORTO SALVO

Tel: +351-214 23 53 50 Fax: +351-214 21 60 36

SOUTH AFRICA

Laboratory Specialities (Pty)

t/a Thermo Fisher Scientific

Thermo Fisher Scientific Building

Unit A, 3 Susan Street
Strijdom Park, Randburg 2169

Tel: +27 11 792 6790 Fax: +27 11 793 1064

SPAIN

Thermo Fisher Diagnostics S.L.U.

Avda. Alcalde Barnils nº 70 Planta 2

Edificio Onada

08174 Sant Cugat del Vallés, Barcelona

Tel: +34-935 765 800 Fax: +34-935 765 820

SWEDEN

Thermo Fisher Diagnostics AB

C/o Phadia AB

P O Box 6460 SE-751 37 UPPSALA

Tel: +46-18 16 60 60 Fax: +46-18 16 63 24

SWITZERLAND

Thermo Fisher Diagnostics AG

Sennweidstrasse 46

CH-6312 STEINHAUSEN

Tel: +41-43 343 40 50 Fax: +41-43 343 40 51

TAIWAN

Phadia Taiwan Inc.

6F-1, No. 85, Jhouzih St., NeiHu District

Taipei, Taiwan, Republic of China

Tel: +886 2 8751 6655 Fax: +886 2 8751 5353

UNITED KINGDOM

Thermo Fisher Diagnostics Ltd

16 Shenley Pavilions, Chalkdell Drive

Shenley Wood, Milton Keynes, MK5 6LB

Tel: +44-1908 76 91 10 Fax: +44-844 324 94 95

USA

Phadia US Inc.

4169 Commercial Avenue

Portage, Michigan 49002

Tel: +1 800-346-4364 (Toll Free) Fax: +1 269 492-7541

OTHER COUNTRIES

Phadia AB,

Distributor Sales

P O Box 6460, SE-751 37 UPPSALA

Tel: +46 18 16 50 00 Fax: +46 18 16 63 65



Phadia AB,
Rapsgatan 7P, P. O. Box 6460, 751 37 Uppsala, Sweden
Tel: +46 18 16 50 00 Fax: +46 18 14 03 58



Issued June 2014. Revised November 2018.
© Phadia AB, Uppsala, Sweden.